

Antibodies Product Datasheet

Product Name:	E4F1 Antibody
Catalog Number:	C43564
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human E4F1
Concentration:	1 mg/ml
Gene Name:	E4F1 E4F
Protein Name:	Transcription factor E4F1 (EC 6.3.2.-) (E4F transcription factor 1) (Putative E3 ubiquitin-protein ligase E4F1) (Transcription factor E4F) (p120E4F) (p50E4F)
Other Name:	
Human Swiss Prot No.:	Q66K89
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human E4F1
Storage and Stability:	-20°C/1 year