

Antibodies Product Datasheet

Product Name:	EBLN2 Antibody
Catalog Number:	C43566
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EBLN2
Concentration:	1 mg/ml
Gene Name:	EBLN2 GK006
Protein Name:	Endogenous Bornavirus-like nucleoprotein 2 (Endogenous Borna-like N element-2) (EBLN-2)
Other Name:	
Human Swiss Prot No.:	Q6P2I7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EBLN2
Storage and Stability:	-20°C/1 year