

Antibodies Product Datasheet

Product Name:	ECHA Antibody
Catalog Number:	C43572
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ECHA
Concentration:	1 mg/ml
Gene Name:	HADHA HADH
Protein Name:	ECHA
Other Name:	Trifunctional enzyme subunit alpha, mitochondrial (78 kDa gastrin-binding protein) (TP-alpha) [Includes: Long-chain enoyl-CoA hydratase (EC 4.2.1.17); Long chain 3-hydroxyacyl-CoA dehydrogenase (EC 1.1.1.211)]
Human Swiss Prot No.:	P40939
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.301% sodium azide.
Specificity:	This antibody detects endogenous levels of ECHA at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year