

Antibodies Product Datasheet

Product Name:	ECHD1 Antibody
Catalog Number:	C43573
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ECHD1
Concentration:	1 mg/ml
Gene Name:	ECHDC1
Protein Name:	Ethylmalonyl-CoA decarboxylase (EC 4.1.1.94) (Enoyl-CoA hydratase domain-containing protein 1) (Methylmalonyl-CoA decarboxylase) (MMCD) (EC 4.1.1.41)
Other Name:	
Human Swiss Prot No.:	Q9NTX5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ECHD1
Storage and Stability:	-20°C/1 year