

Antibodies Product Datasheet

Product Name:	EFS Antibody
Catalog Number:	C43589
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EFS
Concentration:	1 mg/ml
Gene Name:	EFS CASS3
Protein Name:	EFS
Other Name:	Embryonal Fyn-associated substrate (hEFS) (Cas scaffolding protein family member 3)
Human Swiss Prot No.:	O43281
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.86% sodium azide.
Specificity:	This antibody detects endogenous levels of EFS at Human/Mouse
Storage and Stability:	-20°C/1 year