

Antibodies Product Datasheet

Product Name:	EGFL7 Antibody
Catalog Number:	C43590
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	EGFL7 MEGF7 UNQ187/PRO1449
Protein Name:	Epidermal growth factor-like protein 7 (EGF-like protein 7) (Multiple epidermal growth factor-like domains protein 7) (Multiple EGF-like domains protein 7) (NOTCH4-like protein) (Vascular endothelial statin) (VE-statin) (Zneu1)
Other Name:	
Human Swiss Prot No.:	Q9UHF1
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	EGFL7 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year