

Antibodies Product Datasheet

Product Name:	EHD1 Antibody
Catalog Number:	C43593
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EHD1
Concentration:	1 mg/ml
Gene Name:	EHD1 PAST PAST1 CDABP0131
Protein Name:	EH domain-containing protein 1 (PAST homolog 1) (hPAST1) (Testilin)
Other Name:	
Human Swiss Prot No.:	Q9H4M9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EHD1
Storage and Stability:	-20°C/1 year