

Antibodies Product Datasheet

Product Name:	EHMT2 Antibody
Catalog Number:	C43595
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	EHMT2 BAT8 C6orf30 G9A KMT1C NG36
Protein Name:	Histone-lysine N-methyltransferase EHMT2 (EC 2.1.1.-) (EC 2.1.1.43) (Euchromatic histone-lysine N-methyltransferase 2) (HLA-B-associated transcript 8) (Histone H3-K9 methyltransferase 3) (H3-K9-HMTase 3) (Lysine N-methyltransferase 1C) (Protein G9a)
Other Name:	
Human Swiss Prot No.:	Q96KQ7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	EHMT2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year