

Antibodies Product Datasheet

Product Name:	EIF3H Antibody
Catalog Number:	C43600
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EIF3H
Concentration:	1 mg/ml
Gene Name:	EIF3H EIF3S3
Protein Name:	EIF3H
Other Name:	Eukaryotic translation initiation factor 3 subunit H (eIF3h) (Eukaryotic translation initiation factor 3 subunit 3) (eIF-3-gamma) (eIF3 p40 subunit)
Human Swiss Prot No.:	O15372
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.16% sodium azide.
Specificity:	This antibody detects endogenous levels of EIF3H at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year