

Antibodies Product Datasheet

Product Name:	EMC4 Antibody
Catalog Number:	C43614
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EMC4
Concentration:	1 mg/ml
Gene Name:	EMC4 TMEM85 HSPC184 PIG17
Protein Name:	ER membrane protein complex subunit 4 (Cell proliferation-inducing gene 17 protein) (Transmembrane protein 85)
Other Name:	
Human Swiss Prot No.:	Q5J8M3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EMC4
Storage and Stability:	-20°C/1 year