

Antibodies Product Datasheet

Product Name:	EMIL3 Antibody
Catalog Number:	C43616
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EMIL3
Concentration:	1 mg/ml
Gene Name:	EMILIN3 C20orf130 EMILIN5
Protein Name:	EMIL3
Other Name:	EMILIN-3 (EMILIN-5) (Elastin microfibril interface-located protein 3) (Elastin microfibril interfacer 3) (Elastin microfibril interface-located protein 5) (Elastin microfibril interfacer 5)
Human Swiss Prot No.:	Q9NT22
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.313% sodium azide.
Specificity:	This antibody detects endogenous levels of EMIL3 at Human/Mouse
Storage and Stability:	-20°C/1 year