

Antibodies Product Datasheet

Product Name:	ENDO U Antibody
Catalog Number:	C43618
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ENDO U
Concentration:	1 mg/ml
Gene Name:	ENDO U
Protein Name:	Poly(U)-specific endoribonuclease (EC 3.1.-.-) (Placental protein 11) (PP11) (Protein endoU) (Uridylate-specific endoribonuclease)
Other Name:	
Human Swiss Prot No.:	P21128
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ENDO U
Storage and Stability:	-20°C/1 year