

Antibodies Product Datasheet

Product Name:	ENTP2 Antibody
Catalog Number:	C43625
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ENTP2
Concentration:	1 mg/ml
Gene Name:	ENTPD2 CD39L1
Protein Name:	Ectonucleoside triphosphate diphosphohydrolase 2 (NTPDase 2) (EC 3.6.1.-) (CD39 antigen-like 1) (Ecto-ATP diphosphohydrolase 2) (Ecto-ATPDase 2) (Ecto-ATPase 2)
Other Name:	
Human Swiss Prot No.:	Q9Y5L3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ENTP2
Storage and Stability:	-20°C/1 year