

Antibodies Product Datasheet

Product Name:	ERF3B Antibody
Catalog Number:	C43637
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ERF3B
Concentration:	1 mg/ml
Gene Name:	GSPT2 ERF3B
Protein Name:	Eukaryotic peptide chain release factor GTP-binding subunit ERF3B (Eukaryotic peptide chain release factor subunit 3b) (eRF3b) (G1 to S phase transition protein 2 homolog)
Other Name:	
Human Swiss Prot No.:	Q8IYD1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human ERF3B
Storage and Stability:	-20°C/1 year