

Antibodies Product Datasheet

Product Name:	ERR1 Antibody
Catalog Number:	C43641
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 160-240
Concentration:	1 mg/ml
Gene Name:	ESRRA ERR1 ESRL1 NR3B1
Protein Name:	Steroid hormone receptor ERR1 (Estrogen receptor-like 1) (Estrogen-related receptor alpha) (ERR-alpha) (Nuclear receptor subfamily 3 group B member 1)
Other Name:	
Human Swiss Prot No.:	P11474
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ERR1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year