

Antibodies Product Datasheet

Product Name:	ERR2 Antibody
Catalog Number:	C43642
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	ESRRB ERRB2 ESRL2 NR3B2
Protein Name:	Steroid hormone receptor ERR2 (ERR beta-2) (Estrogen receptor-like 2) (Estrogen-related receptor beta) (ERR-beta) (Nuclear receptor subfamily 3 group B member 2)
Other Name:	
Human Swiss Prot No.:	O95718
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ERR2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year