

Antibodies Product Datasheet

Product Name:	ERR3 Antibody
Catalog Number:	C43643
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ESRRG ERR3 ERRG2 KIAA0832 NR3B3
Protein Name:	Estrogen-related receptor gamma (ERR gamma-2) (Estrogen receptor-related protein 3) (Nuclear receptor subfamily 3 group B member 3)
Other Name:	
Human Swiss Prot No.:	P62508
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ERR3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year