

Antibodies Product Datasheet

Product Name:	ESRP2 Antibody
Catalog Number:	C43647
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ESRP2 RBM35B PP7059
Protein Name:	Epithelial splicing regulatory protein 2 (RNA-binding motif protein 35B) (RNA-binding protein 35B)
Other Name:	
Human Swiss Prot No.:	Q9H6T0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ESRP2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year