

Antibodies Product Datasheet

Product Name:	EST1 Antibody
Catalog Number:	C43648
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EST1
Concentration:	1 mg/ml
Gene Name:	CES1 CES2 SES1
Protein Name:	EST1
Other Name:	Liver carboxylesterase 1 (Acyl-coenzyme A:cholesterol acyltransferase) (ACAT) (Brain carboxylesterase hBr1) (Carboxylesterase 1) (CE-1) (hCE-1) (EC 3.1.1.1) (Cocaine carboxylesterase) (Egasyn) (HMSE) (Methylumbelliferyl-acetate deacetylase 1) (EC 3.1.1.56)
Human Swiss Prot No.:	P23141
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.142% sodium azide.
Specificity:	This antibody detects endogenous levels of EST1 at Human/Mouse
Storage and Stability:	-20°C/1 year