

Antibodies Product Datasheet

Product Name:	EST1A Antibody
Catalog Number:	C43649
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EST1A
Concentration:	1 mg/ml
Gene Name:	SMG6 C17orf31 EST1A KIAA0732
Protein Name:	Telomerase-binding protein EST1A (EC 3.1.-.-) (EST1-like protein A) (Ever shorter telomeres 1A) (Smg-6 homolog) (Telomerase subunit EST1A) (hSmg5/7a)
Other Name:	
Human Swiss Prot No.:	Q86US8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EST1A
Storage and Stability:	-20°C/1 year