

Antibodies Product Datasheet

Product Name:	EST3 Antibody
Catalog Number:	C43650
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EST3
Concentration:	1 mg/ml
Gene Name:	CES3 UNQ869/PRO1887
Protein Name:	EST3
Other Name:	Carboxylesterase 3 (EC 3.1.1.1) (Liver carboxylesterase 31 homolog)
Human Swiss Prot No.:	Q6UWW8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.193% sodium azide.
Specificity:	This antibody detects endogenous levels of EST3 at Human
Storage and Stability:	-20°C/1 year