

Antibodies Product Datasheet

Product Name:	EX3L2 Antibody
Catalog Number:	C43661
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EX3L2
Concentration:	1 mg/ml
Gene Name:	EXOC3L2 XTP7
Protein Name:	Exocyst complex component 3-like protein 2 (HBV X-transactivated gene 7 protein) (HBV XAg-transactivated protein 7)
Other Name:	
Human Swiss Prot No.:	Q2M3D2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EX3L2
Storage and Stability:	-20°C/1 year