

Antibodies Product Datasheet

Product Name:	EXTL3 Antibody
Catalog Number:	C43667
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	EXTL3 EXTL1L EXTR1 KIAA0519
Protein Name:	Exostosin-like 3 (EC 2.4.1.223) (EXT-related protein 1) (Glucuronyl-galactosyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase) (Hereditary multiple exostoses gene isolog) (Multiple exostosis-like protein 3) (Putative tumor suppressor protein EXTL3)
Other Name:	
Human Swiss Prot No.:	O43909
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	EXTL3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year