

Antibodies Product Datasheet

Product Name:	EYA2 Antibody
Catalog Number:	C43668
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EYA2
Concentration:	1 mg/ml
Gene Name:	EYA2 EAB1
Protein Name:	Eyes absent homolog 2 (EC 3.1.3.48)
Other Name:	
Human Swiss Prot No.:	O00167
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human EYA2
Storage and Stability:	-20°C/1 year