

Antibodies Product Datasheet

Product Name:	F125A Antibody
Catalog Number:	C43677
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human F125A
Concentration:	1 mg/ml
Gene Name:	FAM125A CFBP
Protein Name:	Multivesicular body subunit 12A (CIN85/CD2AP family-binding protein) (ESCRT-I complex subunit MVB12A) (Protein FAM125A)
Other Name:	
Human Swiss Prot No.:	Q96EY5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human F125A
Storage and Stability:	-20°C/1 year