

Antibodies Product Datasheet

Product Name:	F16A2 Antibody
Catalog Number:	C43689
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human F16A2
Concentration:	1 mg/ml
Gene Name:	FAM160A2 C11orf56 KIAA1759
Protein Name:	F16A2
Other Name:	FTS and Hook-interacting protein (FHIP) (Protein FAM160A2)
Human Swiss Prot No.:	Q8N612
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.288% sodium azide.
Specificity:	This antibody detects endogenous levels of F16A2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year