

Antibodies Product Datasheet

Product Name:	FA156 Antibody
Catalog Number:	C43712
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FA156
Concentration:	1 mg/ml
Gene Name:	FAM156A TMEM29 PP12994 PRO0659; FAM156B TMEM29B
Protein Name:	Protein FAM156A/FAM156B (Transmembrane protein 29)
Other Name:	
Human Swiss Prot No.:	Q8NDB6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FA156
Storage and Stability:	-20°C/1 year