

Antibodies Product Datasheet

Product Name:	FABP5 Antibody
Catalog Number:	C43729
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	FABP5
Protein Name:	Fatty acid-binding protein, epidermal (Epidermal-type fatty acid-binding protein) (E-FABP) (Fatty acid-binding protein 5) (Psoriasis-associated fatty acid-binding protein homolog) (PA-FABP)
Other Name:	
Human Swiss Prot No.:	Q01469
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FABP5 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year