

Antibodies Product Datasheet

Product Name:	FABP6 Antibody
Catalog Number:	C43730
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	FABP6 ILBP ILLBP
Protein Name:	Gastrotropin (GT) (Fatty acid-binding protein 6) (Ileal lipid-binding protein) (ILBP) (Intestinal 15 kDa protein) (I-15P) (Intestinal bile acid-binding protein) (I-BABP)
Other Name:	
Human Swiss Prot No.:	P51161
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FABP6 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year