

Antibodies Product Datasheet

Product Name:	FBSP1 Antibody
Catalog Number:	C43748
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FBSP1
Concentration:	1 mg/ml
Gene Name:	FBXO45 FBX45
Protein Name:	F-box/SPRY domain-containing protein 1 (F-box only protein 45) (hFbxo45)
Other Name:	
Human Swiss Prot No.:	P0C2W1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FBSP1
Storage and Stability:	-20°C/1 year