

Antibodies Product Datasheet

Product Name:	FBX8 Antibody
Catalog Number:	C43754
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FBX8
Concentration:	1 mg/ml
Gene Name:	FBXO8 FBS FBX8 DC10 UNQ1877/PRO4320
Protein Name:	F-box only protein 8 (F-box/SEC7 protein FBS)
Other Name:	
Human Swiss Prot No.:	Q9NRD0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FBX8
Storage and Stability:	-20°C/1 year