

Antibodies Product Datasheet

Product Name:	FCG3A Antibody
Catalog Number:	C43758
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 90-170
Concentration:	1 mg/ml
Gene Name:	FCGR3A CD16A FCG3 FCGR3 IGFR3
Protein Name:	Low affinity immunoglobulin gamma Fc region receptor III-A (CD16a antigen) (Fc-gamma RIII-alpha) (Fc-gamma RIII) (Fc-gamma RIIIa) (FcRIII) (FcRIIIa) (FcR-10) (IgG Fc receptor III-2) (CD antigen CD16a)
Other Name:	
Human Swiss Prot No.:	P08637
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCG3A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year