

Antibodies Product Datasheet

Product Name:	FCN2 Antibody
Catalog Number:	C43762
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	FCN2 FCNL
Protein Name:	Ficolin-2 (37 kDa elastin-binding protein) (Collagen/fibrinogen domain-containing protein 2) (EBP-37) (Ficolin-B) (Ficolin-beta) (Hucolin) (L-ficolin) (Serum lectin p35)
Other Name:	
Human Swiss Prot No.:	Q15485
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCN2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year