

Antibodies Product Datasheet

Product Name:	FCRL1 Antibody
Catalog Number:	C43764
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	FCRL1 FCRH1 IFGP1 IRTA5
Protein Name:	Fc receptor-like protein 1 (FcR-like protein 1) (FcRL1) (Fc receptor homolog 1) (FcRH1) (IFGP family protein 1) (hIFGP1) (Immune receptor translocation-associated protein 5) (CD antigen CD307a)
Other Name:	
Human Swiss Prot No.:	Q96LA6
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCRL1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year