

Antibodies Product Datasheet

Product Name:	FCRL2 Antibody
Catalog Number:	C43765
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	FCRL2 FCRH2 IFGP4 IRTA4 SPAP1 UNQ9236/PRO31998
Protein Name:	Fc receptor-like protein 2 (FcR-like protein 2) (FcRL2) (Fc receptor homolog 2) (FcRH2) (IFGP family protein 4) (Immunoglobulin receptor translocation-associated protein 4) (SH2 domain-containing phosphatase anchor protein 1) (CD antigen CD307b)
Other Name:	
Human Swiss Prot No.:	Q96LA5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCRL2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year