

Antibodies Product Datasheet

Product Name:	FCRL4 Antibody
Catalog Number:	C43766
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	FCRL4 FCRH4 IFGP2 IRTA1
Protein Name:	Fc receptor-like protein 4 (FcR-like protein 4) (FcRL4) (Fc receptor homolog 4) (FcRH4) (IFGP family protein 2) (hIFGP2) (Immune receptor translocation-associated protein 1) (CD antigen CD307d)
Other Name:	
Human Swiss Prot No.:	Q96PJ5
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCRL4 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year