

Antibodies Product Datasheet

Product Name:	FCRLB Antibody
Catalog Number:	C43767
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 340-420
Concentration:	1 mg/ml
Gene Name:	FCRLB FCRL2 FCRLM2 FCRY FREB2
Protein Name:	Fc receptor-like B (Fc receptor homolog expressed in B-cells protein 2) (FREB-2) (Fc receptor-like and mucin-like protein 2) (Fc receptor-like protein 2) (Fc receptor-related protein Y) (FCRY)
Other Name:	
Human Swiss Prot No.:	Q6BAA4
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FCRLB Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year