

Antibodies Product Datasheet

Product Name:	FEZ2 Antibody
Catalog Number:	C43769
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FEZ2
Concentration:	1 mg/ml
Gene Name:	FEZ2
Protein Name:	FEZ2
Other Name:	Fasciculation and elongation protein zeta-2 (Zygin II) (Zygin-2)
Human Swiss Prot No.:	Q9UHY8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.34% sodium azide.
Specificity:	This antibody detects endogenous levels of FEZ2 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year