

Antibodies Product Datasheet

Product Name:	FGD1 Antibody
Catalog Number:	C43770
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGD1
Concentration:	1 mg/ml
Gene Name:	FGD1 FGDY ZFYVE3
Protein Name:	FGD1
Other Name:	FYVE, RhoGEF and PH domain-containing protein 1 (Faciogenital dysplasia 1 protein) (Rho/Rac guanine nucleotide exchange factor FGD1) (Rho/Rac GEF) (Zinc finger FYVE domain-containing protein 3)
Human Swiss Prot No.:	P98174
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.256% sodium azide.
Specificity:	This antibody detects endogenous levels of FGD1 at Human/Mouse
Storage and Stability:	-20°C/1 year