

Antibodies Product Datasheet

Product Name:	FGD4 Antibody
Catalog Number:	C43772
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FGD4
Concentration:	1 mg/ml
Gene Name:	FGD4 FRABP ZFYVE6
Protein Name:	FGD4
Other Name:	FYVE, RhoGEF and PH domain-containing protein 4 (Actin filament-binding protein frabin) (FGD1-related F-actin-binding protein) (Zinc finger FYVE domain-containing protein 6)
Human Swiss Prot No.:	Q96M96
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.315% sodium azide.
Specificity:	This antibody detects endogenous levels of FGD4 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year