

Antibodies Product Datasheet

Product Name:	FHL17 Antibody
Catalog Number:	C43781
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FHL17
Concentration:	1 mg/ml
Gene Name:	FTHL17
Protein Name:	Ferritin heavy polypeptide-like 17 (Cancer/testis antigen 38) (CT38)
Other Name:	
Human Swiss Prot No.:	Q9BXU8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human FHL17
Storage and Stability:	-20°C/1 year