

Antibodies Product Datasheet

Product Name:	FHR1 Antibody
Catalog Number:	C43784
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	CFHR1 CFHL CFHL1 CFHL1P CFHR1P FHR1 HFL1 HFL2
Protein Name:	Complement factor H-related protein 1 (FHR-1) (H factor-like protein 1) (H-factor-like 1) (H36)
Other Name:	
Human Swiss Prot No.:	Q03591
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FHR1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year