

Antibodies Product Datasheet

Product Name:	FHR5 Antibody
Catalog Number:	C43785
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FHR5
Concentration:	1 mg/ml
Gene Name:	CFHR5 CFHL5 FHR5
Protein Name:	FHR5
Other Name:	Complement factor H-related protein 5 (FHR-5)
Human Swiss Prot No.:	Q9BXR6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.165% sodium azide.
Specificity:	This antibody detects endogenous levels of FHR5 at Human
Storage and Stability:	-20°C/1 year