

Antibodies Product Datasheet

Product Name:	FKB1A Antibody
Catalog Number:	C43790
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	FKBP1A FKBP1 FKBP12
Protein Name:	Peptidyl-prolyl cis-trans isomerase FKBP1A (PPIase FKBP1A) (EC 5.2.1.8) (12 kDa FK506-binding protein) (12 kDa FKBP) (FKBP-12) (Calstabin-1) (FK506-binding protein 1A) (FKBP-1A) (Immunophilin FKBP12) (Rotamase)
Other Name:	
Human Swiss Prot No.:	P62942
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FKB1A Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year