

Antibodies Product Datasheet

Product Name:	FKBP9 Antibody
Catalog Number:	C43797
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human FKBP9
Concentration:	1 mg/ml
Gene Name:	FKBP9 FKBP60 FKBP63
Protein Name:	FKBP9
Other Name:	Peptidyl-prolyl cis-trans isomerase FKBP9 (PPIase FKBP9) (EC 5.2.1.8) (63 kDa FK506-binding protein) (63 kDa FKBP) (FKBP-63) (FK506-binding protein 9) (FKBP-9) (Rotamase)
Human Swiss Prot No.:	O95302
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.181% sodium azide.
Specificity:	This antibody detects endogenous levels of FKBP9 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year