

Antibodies Product Datasheet

Product Name:	FLRT3 Antibody
Catalog Number:	C43804
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	FLRT3 KIAA1469 UNQ856/PRO1865
Protein Name:	Leucine-rich repeat transmembrane protein FLRT3 (Fibronectin-like domain-containing leucine-rich transmembrane protein 3)
Other Name:	
Human Swiss Prot No.:	Q9NZU0
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FLRT3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year