

Antibodies Product Datasheet

Product Name:	FOG1 Antibody
Catalog Number:	C43813
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 620-700
Concentration:	1 mg/ml
Gene Name:	ZFPM1 FOG1 ZFN89A
Protein Name:	Zinc finger protein ZFPM1 (Friend of GATA protein 1) (FOG-1) (Friend of GATA 1) (Zinc finger protein 89A) (Zinc finger protein multitype 1)
Other Name:	
Human Swiss Prot No.:	Q8IX07
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FOG1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year