

Antibodies Product Datasheet

Product Name:	FOH1B Antibody
Catalog Number:	C43814
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 370-450
Concentration:	1 mg/ml
Gene Name:	FOLH1B PSMAL GIG26
Protein Name:	Putative N-acetylated-alpha-linked acidic dipeptidase (NAALADase) (EC 3.4.-.-) (Cell growth-inhibiting gene 26 protein) (Prostate-specific membrane antigen-like protein) (Putative folate hydrolase 1B)
Other Name:	
Human Swiss Prot No.:	Q9HBA9
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FOH1B Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year