

Antibodies Product Datasheet

Product Name:	FOXM1 Antibody
Catalog Number:	C43820
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	FOXM1 FKHL16 HFH11 MPP2 WIN
Protein Name:	Forkhead box protein M1 (Forkhead-related protein FKHL16) (Hepatocyte nuclear factor 3 forkhead homolog 11) (HFH-11) (HNF-3/fork-head homolog 11) (M-phase phosphoprotein 2) (MPM-2 reactive phosphoprotein 2) (Transcription factor Trident) (Winged-helix fac
Other Name:	
Human Swiss Prot No.:	Q08050
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FOXM1 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year