

Antibodies Product Datasheet

Product Name:	FPR2 Antibody
Catalog Number:	C43824
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 130-210
Concentration:	1 mg/ml
Gene Name:	FPR2 FPRH1 FPRL1 LXA4R
Protein Name:	N-formyl peptide receptor 2 (FMLP-related receptor I) (FMLP-R-I) (Formyl peptide receptor-like 1) (HM63) (Lipoxin A4 receptor) (LXA4 receptor) (RFP)
Other Name:	
Human Swiss Prot No.:	P25090
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	FPR2 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year